


RespiClick® Device

Patient Teaching Guide


1

Flip open the cap of the RespiClick® inhaler until you hear a click. The dose is now loaded and ready to inhale. Be careful to keep the device upright until after inhalation.




2

Exhale completely and fully away from the inhaler.



3

Close your lips around the mouth-piece, being careful to keep the inhaler upright and not to cover the air vent above the mouthpiece. Inhale sharply and deeply. Hold your breath for as long as you can, up to 10 seconds.



4

Flip cap closed until you hear a click.

